

# MERCK

## Microlab 200

### DIMENSIONS AND WEIGHT

Width	380 mm
Depth	390 mm
Height	140 mm
Weight approx.	7 kg

### INSTRUMENT LOCATION

The Merck Microlab 200 is intended for indoor use. It should be located on a flat surface, which is not exposed to vibration (i.e. centrifuges, etc.) and direct sunlight.

### POWER REQUIREMENTS

Line Voltage	110..240 V nominal, tolerance 10%
Line Frequency	50/60 Hz
Power Consumption	100 VA
Installation category	II

The Merck Microlab 200 is, without any switching, suitable for the whole specified voltage range of 110...240 V nominal.

### ENVIRONMENTAL REQUIREMENTS

- Ambient temperature 10 to 35 °C
- Maximum relative humidity 80% at temperatures up to 31 °C, linear decreasing to 65% at 35 °C
- Altitude up to 3000 m
- Pollution degree 2 (in accordance with IEC 664)



Tel. 5558199705

excelab2014@hotmail.com